

**L13****Simultaneous measurement of sarcomere shortening and calcium transients in primary rat cardiomyocytes exposed to electrical pulses*****Vid Jan, Marko Stručić****University of Ljubljana, Faculty of Electrical Engineering*

Duration of the experiments: 90 min

Max. number of participants: 3

Location: Cell Culture Laboratory 1

Level: Advanced

**PREREQUISITES**

Participants should be familiar with Laboratory safety (S1). Basic knowledge of handling with cells is required for this laboratory practice.

**The aim** of this laboratory practice is to observe the effects of electroporation – using monopolar single electric pulse or a burst of bipolar high frequency irreversible electroporation pulses (H-FIRE) on  $\text{Ca}^{2+}$  homeostasis, contractility and recovery of adult rat ventricular cardiomyocytes.


**THEORETICAL BACKGROUND**

Despite a long history of using electric pulses on the cardiac tissue for defibrillation and pacing, effects of high-voltage (supraphysiological) electric pulses at the level of a single cardiomyocyte (or even neurons and excitable cells in general, for that matter), remain relatively unexplored. It is known that during electroporation, cell membrane is damaged, its conductivity is increased, and cell homeostasis disrupted [1]. Voltage-gated channels can also be affected, although the consensus on the mechanisms and the consequences of electroporation on voltage-gated channels is not yet reached. Interestingly, in a recent molecular dynamics study, the voltage sensor area of the voltage-gated channels was shown to be permanently damaged under the influence of high electric fields [2]. During electroporation,  $\text{Ca}^{2+}$  enters the cell in an uncontrolled way through the permeabilized membrane. Thus, the precise calcium balance, needed for normal cardiomyocyte function is disrupted [3]. This can lead to undesired and unpredictable behaviours. For example, excitation-contraction coupling (pairing cardiomyocyte depolarization with mechanical contraction) depends largely on  $\text{Ca}^{2+}$  signalling and disturbances of  $\text{Ca}^{2+}$  homeostasis can severely impact muscle physiology and give rise to various myopathies and cardiac disorders. Under physiological conditions, cardiomyocytes can efficiently remove the excess intracellular  $\text{Ca}^{2+}$ , which enters during the action potential plateau phase to prepare for the next cycle. However, with electroporation, the uptake of  $\text{Ca}^{2+}$  can be much higher than in normal working conditions; moreover, the cell membrane is damaged and does not reseal for minutes after the treatment. Furthermore, ATP, necessary for membrane resealing and operation of pumps, following electroporation leaks out of the cell [4].

The efficiency of electroporation depends on the parameters of the electric field, cell size, geometry, position, and physiological characteristics of the medium surrounding the cell. Different fluorescent dyes (e.g. Rhod-2, Fluo-4, Fura-2, Fura-3,...) can be used to investigate the influence of these parameters on electroporation and the same dyes can be used to monitor how electroporation affects Ca<sup>2+</sup> transients [5].

## EXPERIMENT

We will simultaneously monitor sarcomere shortenings, i.e., contraction and Ca<sup>2+</sup> transients of rat cardiomyocytes with IonOptix-upgraded Zeiss Axiovert microscope which applies Fast Fourier Transform to measure sarcomere length while at the same time measuring signal from fluorescent Ca<sup>2+</sup> sensitive indicator Rhod-2 AM. As cardiomyocytes are paced, Ca<sup>2+</sup> concentrations inside the cells change, due to its importance in action potentials, while supraphysiological electric pulses, i.e. electroporation, cause changes in Ca<sup>2+</sup> concentrations due to perturbations in membrane permeability (due to electroporation). Once Ca<sup>2+</sup> enters the sarcoplasm of cardiomyocytes, it binds to the dye and changes its fluorescence intensity more than 100-fold (Figure 1).




**Figure 1:** (A) Fluorescence excitation (blue) and emission (red) spectra of Rhod-2 AM for different concentrations of Ca<sup>2+</sup> (image from <https://www.thermofisher.com/order/fluorescence-spectraviewer?SID=srch-svtool&UID=14220ca#!/>). (B) Schematic of the experimental setup (adapted from Chaigne, S. et al., 2022 [3]). (C) Simultaneous acquisition of calcium transients and sarcomere shortenings with IonOptix system.

**Protocol:** The experiments will be performed on adult rat ventricular cardiomyocytes, isolated from female Wistar rats, weighing 250-350 grams [6,7]. Before measurements, cardiomyocytes will be incubated in Tyrode's buffer with 2  $\mu$ M Rhod-2, AM for 20 min at 37 °C (Thermo Fisher Scientific). After the incubation time, cardiomyocytes will be washed with fresh Tyrode's to remove excess dye and then seeded on 25 mm square coverslips for imaging.

Insert fibronectin/gelatin-coated glass coverslips into IonOptix stimulation chamber on fluorescence microscope (Zeiss Axiovert 200) and use 63 $\times$ LD objective. Wait 5 min and then start the perfusion of Tyrode's buffer at 37 °C. IonOptix MyoCam-S3™ Fast CMOS camera will be used to capture

sarcomere shortening while IonOptix photomultiplier tube will be used to capture signals for Rhod-2AM. Measurements will be recorded and processed using IonOptix IonWizard software.

Using an IonOptix MyoPacer Field Stimulator pace cardiomyocytes on the coverslip with bipolar pulses at 2 Hz and 2 ms (2 ms positive phase followed immediately by 2 ms negative phase). IonOptix stimulation chamber is equipped with 2 platinum electrodes (4 mm spacing edge to edge) and our setup also includes a switch box to alternate between pulses from the IonOptix MyoPacer Field Stimulator for physiological pacing and HF-IRE pulse generator (University of Ljubljana) for monopolar or bipolar supraphysiological pulses. Using a framing adapter, frame a cardiomyocyte that contracts in sync with delivered pacing pulses. Measure basal changes in sarcomere length and  $\text{Ca}^{2+}$  transients (capture at least 10 contractions). Afterwards use HF-IRE pulse generator and deliver one monopolar electric pulse of 100  $\mu\text{s}$  or a burst of H-FIRE pulses (25 bipolar pulses, 2  $\mu\text{s}$  positive and negative phase, 2  $\mu\text{s}$  interpulse and interphase delay), which will represent an electroporation (EP) pulse. Start with 100 V/cm for monopolar pulses, and with 300 V/cm for H-FIRE pulses. 1 minute after delivered pulse, start pacing the cell again with 2 ms bipolar pulses. If cell still contracts, deliver another pulse with higher electric field (with increments of 50 V/cm). Continue in the same fashion until cell hypercontracts and does not respond to pacing 1 min after EP pulse delivery. With this approach we expect to observe different types of responses to EP protocol at different electric field strengths (shown in Figure 2). Compare results for sarcomere shortening and  $\text{Ca}^{2+}$  transients of cells treated with monopolar or H-FIRE pulses.



**Figure 2:** (A-D) Simultaneous measurements of calcium transients ( $[\text{Ca}^{2+}]_i$ ), and sarcomere shortening. Primary rat cardiomyocyte was firstly paced at 2 Hz with 8 V bipolar pulses of 2 ms (Pacing). After that the cell will be exposed to either one 100  $\mu\text{s}$  monopolar pulse or one burst of bipolar H-FIRE pulses (25 bipolar pulses, 2  $\mu\text{s}$  positive and negative phase, 2  $\mu\text{s}$  interpulse and interphase delay) (EP protocol). Each pulse or burst of pulses is marked with an asterisk (EP protocol). 1 min after EP protocols 2 ms pacing continued (Pacing after EP protocol). Normally we see different types of responses to electroporation, which correlate with the intensity of electric field at which EP protocol is performed. (A) At lowest electric fields there is no response to EP protocols. (B) When electric field strength is increased EP protocol

induces a contraction of a treated cardiomyocyte. (C) When electric field strength is increased even further an EP pulse induces several spontaneous contractions and oscillations in calcium, but the cell can still be paced 1 min after EP pulse delivery. (D) At highest electric fields spontaneous contractions and high levels of calcium are observed and the cell becomes hypercontracted after EP pulse delivery. Cell is unresponsive to pacing 1 min after EP pulse delivery.

**REFERENCES:**

1. Tung L., Tovar O., Neunlist M., Jain S.K., O'Neill R.J. Effects of strong electrical shock on cardiac muscle tissue. *Ann NY Acad Sci*, 720:160-175, 1994.
2. Rems L., Kasimova M.A., Testa I., Delemotte L. pulsed electric fields can create pores in the voltage sensors of voltage-gated ion channels. *Biophys J*, 119:190-205, 2020.
3. Chaigne S., Sigg D.C., Stewart M.T., Hocini M., Batista Napotnik T., Miklavčič D., Bernus O., Benoist D. Reversible and irreversible effects of electroporation on contractility and calcium homeostasis in isolated cardiac ventricular myocytes. *Circ Arrhythm Electrophysiol*, 15: e011131, 2022.
4. Polajzer T., Jarm T., Miklavčič D. Analysis of damage-associated molecular pattern molecules due to electroporation of cells in vitro. *Radiol Oncol*, 54:317-328, 2020.
5. Batista Napotnik T., Miklavčič, D. In vitro electroporation detection methods – An overview. *Bioelectrochemistry*, 120:166-182, 2018.
6. Louch W.E., Sheehan K.A., Wolska B.M. Methods in cardiomyocyte isolation, culture, and gene transfer. *J Mol Cell Cardiol*, 51:288-298, 2011.
7. Nippert F., Schreckenberger R., Schlüter K.D. Isolation and cultivation of adult rat cardiomyocytes. *J Vis Exp*, 56634, 2017.

**NOTES & RESULTS**

---

