

L6 *E. coli* inactivation by pulsed electric fields in a continuous mode

Saša Haberl Meglič, Karel Flisar

University of Ljubljana, Faculty of Electrical Engineering

Duration of the experiment: day 1: 90 min; day 2: 30 min

Max. number of participants: 4

Location: Microbiological laboratory

Level: Basic

PREREQUISITES

Participants should be familiar with Laboratory safety (S1) and Electroporation hardware safety (S2). No other specific knowledge is required for this laboratory practice.

The aim of this laboratory practice is to inactivate bacteria with electroporation in a flow-through system.

THEORETICAL BACKGROUND

Electroporation (under this name) was first described 50 years ago [1] and causes transient increase in permeability of the cell membrane by applying high-voltage electric field pulses. At stronger electric fields, cells are damaged, leading to cell death. Such application also known as irreversible electroporation has been previously used to inactivate bacteria in water environment [2]. The method gained ground as a tool for microbial inactivation and the influence of different electroporation parameters on microbial viability was extensively studied on various microorganisms [3].

Since PEF microbial inactivation in controlled laboratory conditions showed promise, the idea arose of also removing pathogenic microorganisms from various water sources, hospital wastewaters and liquid food, without destroying vitamins or affecting the food's flavor, color or texture [4-6]. To enable electroporation on a large scale, the development of flow-through processes has been proposed [7]. Thus, a flow-through treatment system consists of a pulse generator that provides continuous pulse treatment, flow-through chambers with electrodes, and a fluid handling system.

Several parameters have been described, which can influence inactivation of microbial cells. Specifically in a continuous flow system the flow rate of a liquid must be adjusted for each bacterial cell to be exposed to appropriate pulse conditions [8].

EXPERIMENT

We will inactivate *Escherichia coli* K12 TOP10 cells carrying plasmid pEGFP-N1, which encodes kanamycin resistance (Clontech Laboratories Inc., Mountain View, CA, USA) in a continuous flow system (see Figure 1) using various electrical pulse parameters.

A prototype square wave pulse generator will be used to generate electrical pulses. The pulses will be monitored using an oscilloscope (LeCroy 9310C). The degree of inactivation will be determined using the plate counting method.




Figure 1: Continuous mode electroporation system. The circuit system includes a flow chamber with electrodes and a prototype square wave pulse generator. Voltage and current are both monitored throughout the experiment.

Protocol 1/2 (Electroporation of bacteria in a continuous mode): Bacterial cells will be grown prior experiment for 17 hours at 37°C in Luria Broth (LB) medium (Sigma-Aldrich Chemie GmbH, Deisenhofen, Germany) with shaking. *E. coli* cells will be then centrifuged (4248 g, 30 min, 4°C) and the pellet will be resuspended in 250 mM sucrose.

The exposure of cells to electrical pulses in a flow chamber in a continuous mode system depends on the geometry of the chamber and the frequency of the pulses at which the electroporator operates. The number of pulses is determined by equation 1. At this flow rate, the desired number of pulses is applied to the fluid and thus to the cells in the flow chamber. Since the volume of our cross-field chamber between the electrodes and the frequency are constant, the flow through the chamber can be determined:

$$q = \frac{f}{n} \cdot Q \quad (1)$$

where q (L/min) is the flow rate, Q (L) is the volume between the two electrodes, and n is the number of pulses received by the fluid in the chamber during the dwell time. For a frequency of 10 Hz, calculate the flow rate (q) at which all the fluid is exposed to at least one pulse. The bacterial cells will be pumped through the system at the calculated flow rate, and pulses will be applied by the prototype pulse generator.

After electroporation treatment, collect 20 µl of the treated sample and prepare dilutions ranging from 10^{-1} to 10^{-6} . Pipette three 10-µl drops of the different dilutions onto LB agar.

To determine the number of bacterial cells in our sample, make serial dilutions of the (untreated) bacterial sample ranging from 10^{-1} to 10^{-7} (dilute 20 µl of the untreated bacterial sample in 180 µl of 0.9% NaCl). Pipette three 10-µl drops of dilutions 10^{-5} , 10^{-6} and 10^{-7} onto LB agar.

Protocol 2/2 (Determining bacterial viability): After 24 hours of incubation at 37°C, count the colony forming units. Viability is expressed as $\log(N/N_0)$, where N is the number of colony forming units per ml in a treated sample and N_0 is the number of colony forming units per ml in an untreated sample.

Example of determining bacterial viability:

You counted 20 CFU in a control sample (dilution 10^{-7}) and 10 CFU in a treated sample (dilution 10^{-5}).

Number of bacterial cells per ml (control sample) = 20×10^7 (dilution factor of sample) $\times 100$ (dilution factor of plating) = 2×10^{10} bacterial cells/ml

Number of bacterial cells per ml (treated sample) = 10×10^5 (dilution factor of sample) $\times 100$ (dilution factor of plating) = 1×10^8 bacterial cells/ml

$\log N/N_0 = \log(1 \times 10^8 / 2 \times 10^{10}) = -2.301$

REFERENCES:

- [1] Neumann E., Rosenheck K. Permeability changes induced by electric impulses in vesicular membranes. *J Membr Biol*, 10:279-90, 1972.
- [2] Gusbeth C., Frey W., Volkmann H., Schwartz T., Bluhm H. Pulsed electric field treatment for bacteria reduction and its impact on hospital wastewater. *Chemosphere*, 75: 228-233, 2009.
- [3] Zhou J.F., Hung Y.C., Xie X. Application of electric field treatment (EFT) for microbial control in water and liquid food. *J Hazard Mater*, 445: 130561, 2023.
- [4] Lytras F., Psakis G., Gatt R., Cebrián G., Raso J., Valdramidis V. Exploring the efficacy of pulsed electric fields (PEF) in microbial inactivation during food processing: A deep dive into the microbial cellular and molecular mechanisms. *IFSET*, 95: 103732, 2024.
- [5] Gomez B., Munekata P.E.S., Gavahian M., Barba F.J., Marti-Quijal F.J., Bolumar T., *et al.* Application of pulsed electric fields in meat and fish processing industries: An overview. *Food Res Int*, 123:95-105, 2019.
- [6] Zhou J., Hung Y-C, Xie X. Application of electric field treatment (EFT) for microbial control in water and liquid food. *J Hazard Mater*, 445:130561, 2023.
- [7] Flisar K., Haberl Meglic S., Morelj J., Golob J., Miklavčič D. Testing a prototype pulse generator for a continuous flow system and its use for *E. coli* inactivation and microalgae lipid extraction. *Bioelectrochemistry*, 100:44-51, 2014.
- [8] Pataro G., Senatore B., Donsi G., Ferrari G. Effect of electric and flow parameters on PEF treatment efficiency. *J Food Eng*, 105:79-88, 2011.

EXPECTED RESULTS




Figure 2: Expected reduction of *E. coli* in a continuous mode electroporation system. Bacterial cells were exposed to 10 or 30 rectangular pulses with a pulse width of 100 µs or 1 ms, a pulse frequency of 10 Hz, and a pulse amplitude of 1500 V, corresponding to a voltage-to-distance ratio of 0.6 kV/cm.

NOTES & RESULTS

